

HOW TO INSERT BONE.TXT

1. click on START
2. click on PROGRAMS
3. click on IMCSRS V 3.2 for Win 95 NT
4. click on imcsdata
5. enter 'imcs' for name, enter 'imcs' for password
6. click on Administrative Menu
7. click on Update Reportable Equipment
8. Put floppy disk containing MMDF in 'a' drive
9. change 'c:\' to read 'a:\'
10. click on okay
11. click on 'bone.txt' to highlight
12. click on okay
13. click on okay (wait)
14. click on ok
15. click ok (transfer complete)
16. click ok (Reportable equipment file was upgraded)
17. click on EXIT (IMCSRS)
18. click on "you will Exit IMCSRS"